Steripack Instrument Case and Tray Systems

INTENDED USE	<ul style="list-style-type: none"> Steripack Case and Tray systems are utilized to secure medical instruments during transport, storage and sterilization processing. 										
INTENDED USER PROFILE	<ul style="list-style-type: none"> Hospital and surgical staff having adequate training and familiarity with the handling of instruments including but not limited to, the loading and unloading of, transport, storage, and sterilization of instruments in a case and tray system prior to and following surgical procedures. 										
DEVICE DESCRIPTION	<ul style="list-style-type: none"> Steripack sterilization systems are manufactured using stainless steel, anodized aluminum, nylon 11, and medical grade silicone materials. Cases and trays are supplied NON-STERILE and must be inspected, cleaned and sterilized before use. Devices are not implantable. 										
WARNINGS 	<ul style="list-style-type: none"> Devices should be reprocessed as soon as possible following use. Instruments must be cleaned separately from cases and trays. If a device is/was used in a patient with, or suspected of having Creutzfeldt-Jakob Disease (CJD), the device cannot be reused and must be destroyed due to the inability to reprocess or sterilize to eliminate the risk of cross-contamination. 										
CAUTION  ONLY	<p>Federal U.S. Law restricts this device to sale, distribution, and use, by, or on order of a physician.</p>										
LIMITATIONS ON REPROCESSING	<p>Repeated processing per this instruction for use has minimal effect on these instruments. End of life is normally determined by wear and damage due to use.</p>										
DISCLAIMER	<p>It is the responsibility of the reprocessor to ensure reprocessing is performed using equipment, materials and personnel in the reprocessing facility and achieves the desired result. This requires validation and routine monitoring of the process. Any deviation by the reprocessor from the instructions provided must be properly evaluated for effectiveness and potential adverse consequences.</p>										
MANUAL CLEANING	<ol style="list-style-type: none"> After each use, wash the sterilization tray with an aluminum safe, neutral pH detergent to avoid faded surface colors and deterioration of the anodized surface. Thoroughly rinse the sterilization tray with warm tap water for 1 minute and dry with a soft, absorbent cloth. 										
INSPECTION	<ul style="list-style-type: none"> Visually inspect devices for damage or wear. Do not use a sterilization tray if the lid does not attach to the base securely or if the sterilization tray is or appears damaged. 										
PACKAGING	<ul style="list-style-type: none"> To maintain sterility, the sterilization tray must be wrapped in an approved sterilization wrap using the AAMI double wrap method or equivalent prior to sterilization. The end user should consult ANSI/AAMI ST79 for additional information on steam sterilization. 										
STERILIZATION	<p>Sterilize with steam. The following are minimum cycles required for steam sterilization of the case assembly with up to 896 grams of metal instruments:</p> <p>Double Wrapped Instrument Case:</p> <table border="1" data-bbox="396 1642 1305 1705"> <thead> <tr> <th>Cycle Type</th> <th>Temperature</th> <th>Exposure Time</th> <th>Pulses</th> <th>Drying Time</th> </tr> </thead> <tbody> <tr> <td>Prevacuum</td> <td>132°C (270°F)</td> <td>4 minutes</td> <td>3</td> <td>20 minutes</td> </tr> </tbody> </table> <p>NOTE: The cases have been validated at these parameters with an 896 gram load of metal instruments. Deviations from this case loading or instrument type will require re-validation by the user.</p> <ul style="list-style-type: none"> Time and temperature parameters required for sterilization vary according to type of sterilizer, cycle design, and packaging material. It is critical that process parameters be validated for each facility's individual type of sterilization equipment and product load configuration. 	Cycle Type	Temperature	Exposure Time	Pulses	Drying Time	Prevacuum	132°C (270°F)	4 minutes	3	20 minutes
Cycle Type	Temperature	Exposure Time	Pulses	Drying Time							
Prevacuum	132°C (270°F)	4 minutes	3	20 minutes							

	<ul style="list-style-type: none"> Only steam sterilization cycles with the above listed parameters have been validated for use and have been shown to be compatible with the device design and specified instrument loading. A facility may choose to use different steam sterilization cycles other than the cycle suggested if the facility has properly validated the cycle to ensure adequate steam penetration and contact with the devices for sterilization.
STORAGE	<ul style="list-style-type: none"> After sterilization, instrument cases should remain in sterilization packaging and be stored in a clean, dry cabinet or storage case. Care should be taken when handling devices to avoid damaging the sterile barrier.
CONTACT	<div style="display: flex; justify-content: space-between; align-items: flex-start;"> <div style="width: 45%;">  <p>Manufactured by: Avalign Technologies 8727 Clinton Park Drive Fort Wayne, IN 46825 1-877-289-1096 www.avalign.com product.questions@avalign.com</p> </div> <div style="width: 15%; text-align: center;"> <div style="border: 1px solid black; padding: 2px; display: inline-block;">EC</div> <div style="border: 1px solid black; padding: 2px; display: inline-block; margin-left: 5px;">REP</div> </div> <div style="width: 35%; text-align: right;">  <p>Authorized Representative: mdi Europa GmbH Langenhagener Str. 71 D-30855 Langenhagen +49-511-3908 9530 www.mdi-europa.com</p> </div> </div>

Label Glossary

Symbol	Title and Translations
	Manufacturer
	Authorized Representative in the European Community
	Lot Number / Batch Code
	Catalogue Number
	Consult Instructions for Use
	Caution
	Federal Law (USA) restricts this device to sale by or on the order of a physician

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