Steripack Instrument Case and Tray Systems

[Translations available](#) (click on language to advance to translation):

LANGUAGES
English (EN)
French (Français - FR)
German (Deutsch - DE)
Italian (Italiano - IT)
Portuguese (Português – PT)
Spanish (Español – ES)
Danish (Dansk – DA)
Dutch (Nederlands - NL)
Greek (Ελληνικά - EL)
Japanese (日本語 – JP)
Swedish (Svenska – SE)

Steripack Instrument Case and Tray Systems

INTENDED USE	<ul style="list-style-type: none">Steripack Case and Tray systems are utilized to secure medical instruments during transport, storage and sterilization processing.															
INTENDED USER PROFILE	<ul style="list-style-type: none">Hospital and surgical staff having adequate training and familiarity with the handling of instruments including but not limited to, the loading and unloading of, transport, storage, and sterilization of instruments in a case and tray system prior to and following surgical procedures.															
INDICATIONS	<ul style="list-style-type: none">The Steripack case and tray system are intended for use in healthcare facilities to organize, enclose, sterilize, transport, and store medical devices and other instrumentation between surgical and other medical uses.The Steripack case and tray system are not intended on their own to maintain sterility; they are intended to be used in conjunction with a legally marketed, validated, FDA=cleared sterilization wrap.Sterilization validation for the worst-case Steripack cases and tray system included surgical instrument such as rongeurs, forceps, wrenches, cutters, pliers, etc. The Steripack cases and tray system were validated for with up to a 9.35 lb (4.24 kg) load of metal instruments and polymer handled instruments. <p>Sterilization Parameters:</p> <table><tr><td>Cycle Type</td><td>Temperature</td><td>Exposure Time</td><td>Pulses</td><td>Drying Time</td></tr><tr><td>Prevacuum</td><td>132°C (270°F)</td><td>4 minutes</td><td>3</td><td>40 minutes</td></tr><tr><td>Prevacuum</td><td>134°C (273°F)</td><td>3 minutes</td><td>3</td><td>45 minutes</td></tr></table>	Cycle Type	Temperature	Exposure Time	Pulses	Drying Time	Prevacuum	132°C (270°F)	4 minutes	3	40 minutes	Prevacuum	134°C (273°F)	3 minutes	3	45 minutes
Cycle Type	Temperature	Exposure Time	Pulses	Drying Time												
Prevacuum	132°C (270°F)	4 minutes	3	40 minutes												
Prevacuum	134°C (273°F)	3 minutes	3	45 minutes												
DEVICE DESCRIPTION	<ul style="list-style-type: none">Steripack sterilization systems are manufactured using stainless steel, anodized aluminum, nylon 11, and medical grade silicone materials.Cases and trays are supplied NON-STERILE and must be inspected, cleaned and sterilized before use.Devices are not implantable.															
WARNINGS <div></div>	<ul style="list-style-type: none">Devices should be reprocessed as soon as possible following use. Instruments must be cleaned separately from cases and trays.If a device is/was used in a patient with, or suspected of having Creutzfeldt-Jakob Disease (CJD), the device cannot be reused and must be destroyed due to the inability to reprocess or sterilize to eliminate the risk of cross-contamination.															
CAUTION <div></div> <div>ONLY</div>	Federal U.S. Law restricts this device to sale, distribution, and use, by, or on order of a physician.															
LIMITATIONS ON REPROCESSING	Repeated processing per this instruction for use has minimal effect on these instruments. End of life is normally determined by wear and damage due to use.															
DISCLAIMER	It is the responsibility of the reprocessor to ensure reprocessing is performed using equipment, materials and personnel in the reprocessing facility and achieves the desired result. This requires validation and routine monitoring of the process. Any deviation by the reprocessor from the instructions provided must be properly evaluated for effectiveness and potential adverse consequences.															
MANUAL CLEANING	<ol style="list-style-type: none">After each use, wash the sterilization tray with a soft sponge and an aluminum safe, neutral pH detergent. A neutral pH detergent is required to avoid faded surface colors and deterioration of the anodized surface.Thoroughly rinse the sterilization tray with warm tap water for 1 minute and dry with a soft, absorbent cloth.															
INSPECTION	<ul style="list-style-type: none">Visually inspect devices for damage or wear.Do not use a sterilization tray if the lid does not attach to the base securely or if the sterilization tray is or appears damaged.															
PACKAGING	<ul style="list-style-type: none">To maintain sterility, the sterilization tray must be wrapped in an approved sterilization wrap using the AAMI double wrap method or equivalent prior to sterilization.The end user should consult ANSI/AAMI ST79 for additional information on steam sterilization.															

STERILIZATION	<p>Sterilize with steam. The following are minimum cycles required for steam sterilization of the case assembly with up to a 9.35 lb (4.24 kg) load of metal instruments and polymer handled instruments:</p> <p>Double Wrapped Instrument Case:</p> <table><tr><th>Cycle Type</th><th>Temperature</th><th>Exposure Time</th><th>Pulses</th><th>Drying Time</th></tr><tr><td>Prevacuum</td><td>132°C (270°F)</td><td>4 minutes</td><td>3</td><td>40 minutes</td></tr><tr><td>Prevacuum</td><td>134°C (273°F)</td><td>3 minutes</td><td>3</td><td>45 minutes</td></tr></table> <p>NOTE: The cases have been validated at these parameters with an 9.35 lb (4.24 kg) load of metal instruments and polymer handled instruments. Deviations from this case loading or instrument type will require re-validation by the user.</p> <ul style="list-style-type: none">Time and temperature parameters required for sterilization vary according to type of sterilizer, cycle design, and packaging material. It is critical that process parameters be validated for each facility’s individual type of sterilization equipment and product load configuration.Only steam sterilization cycles with the above listed parameters have been validated for use and have been shown to be compatible with the device design and specified instrument loading. A facility may choose to use different steam sterilization cycles other than the cycle suggested if the facility has properly validated the cycle to ensure adequate steam penetration and contact with the devices for sterilization.Water droplets and visible signs of moisture on sterile packaging or the tape used to secure it, may compromise sterility of processed loads or be indicative of a sterilization process failure. Visually check outside wrapper for dryness. If there are water droplets or visible moisture on the exterior of the package or on the tape used to secure it, the pack or instrument tray is considered unacceptable. Repackage and re-sterilize sterilization packages with visible signs of moisture.	Cycle Type	Temperature	Exposure Time	Pulses	Drying Time	Prevacuum	132°C (270°F)	4 minutes	3	40 minutes	Prevacuum	134°C (273°F)	3 minutes	3	45 minutes
Cycle Type	Temperature	Exposure Time	Pulses	Drying Time												
Prevacuum	132°C (270°F)	4 minutes	3	40 minutes												
Prevacuum	134°C (273°F)	3 minutes	3	45 minutes												
STORAGE	<ul style="list-style-type: none">After sterilization, instrument cases should remain in sterilization packaging and be stored in a clean, dry cabinet or storage case.Care should be taken when handling devices to avoid damaging the sterile barrier.															
CONTACT	<p>Notice to Patient and User: Any serious incident that has occurred in relation to the medical devices should be reported to the manufacturer and the competent authority of the EU Member State in which the user and/or patient is established.</p> <div><div><p>Manufactured by: Avalign Technologies 8727 Clinton Park Drive Fort Wayne, IN 46825 1-877-289-1096 www.avalign.com product.questions@avalign.com</p></div><div><div><div>EC</div><div>REP</div></div><p>Authorized Representative:</p><p>Emergo Europe, Prinsessegracht 20, 2514 AP The Hague, The Netherlands</p></div><div></div></div>															

Label Glossary

Symbol	Title and Translations	Symbol	Title and Translations
	Manufacturer and Date of Manufacture		Consult Instructions for Use
	Authorized Representative in the European Community		Caution
	Lot Number / Batch Code		Federal Law (USA) restricts this device to sale by or on the order of a physician
	Catalogue Number		Medical Device

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