## Tissue and Palatal Retractor: Instrumentation Instructions for Use

Translations available (click on language to advance to translation):

| LANGUAGE                    |
|-----------------------------|
| English (EN)                |
| French (Français - FR)      |
| German (Deutsch - DE)       |
| Italian (Italiano - IT)     |
| Portuguese (Português – PT) |
| Spanish (Español – ES)      |
| Bulgarian (Български – BG)  |
| Croatian (Hrvatski -HR)     |
| Czech (Česky – CS)          |
| Danish (Dansk – DA)         |
| Dutch (Nederlands - NL)     |
| Estonian (Eesti – ET)       |

| LANGUAGE                      |
|-------------------------------|
| Finnish (Finnish - FIN)       |
| Greek (Ελληνικά - EL)         |
| Hungarian (Magyar – HU)       |
| Latvian (Latviski – LV)       |
| Lithuanian (Lietuvių k. – LT) |
| Norwegian (Norsk – NO)        |
| Polish (Polski – PL)          |
| Romanian (Română – RO)        |
| Slovak (Slovenčina – SK)      |
| Slovenian (Slovenščina – SL)  |
| Swedish (Svenska – SE)        |
| Turkish (Türkçe – TR)         |

### Part List

| Catalog Number | Description   | UDI-DI         |
|----------------|---|----------------|
| 03.520.000     | Double Ended Retractor                                | 00190776028138 |
| 03.520.001     | Single Hook Skin Retractor, 3mm                       | 00190776028145 |
| 03.520.002     | Double Hook Skin Retractor, 2mm                       | 00190776028152 |
| 03.520.004     | Palatal Retractor Frame                               | 00190776121167 |
| 03.520.005     | Palatal Retractor Frame Blade Short                   | 00190776121174 |
| 03.520.006     | Palatal Retractor Frame Blade Medium                  | 00190776121181 |
| 03.520.007     | Palatal Retractor Frame Blade Long                    | 00190776121198 |
| 03.520.008     | Double Ended Periosteal Elevator, Straight, 19 cm     | 00190776028169 |
| 03.520.009     | Double Ended Periosteal Elevator, Angle, 19 cm        | 00190776028176 |
| 03.520.010     | Long Tissue Forceps, DeBakey Style                    | 00190776028183 |
| 03.520.011     | Tissue Forceps, Adson Style                           | 00190776028190 |
| 03.520.012     | Tissue Forceps, Adson Brown Style                     | 00190776028206 |
| 03.520.013     | Scissors, Kelly Style                                 | 00190776028213 |
| 03.520.015     | Scissors, Gorney Style, Long, Curved                  | 00190776028237 |
| 03.520.016     | Scissors, Flat Pointed, Curved                        | 00190776028244 |
| 03.520.017     | Scissors, Iris Style, Curved                          | 00190776028251 |
| 03.520.018     | Needle Holders, Tip Crile-Wood Style, Short, Smooth   | 00190776028268 |
| 03.520.019     | Needle Holder, Crile-Wood Style, Medium               | 00190776028275 |
| 03.520.020     | Needle Holder, Necked, Long                           | 00190776028282 |
| 03.520.021     | Osteotome, Epker Style, 4mm                           | 00190776028299 |
| 03.520.022     | Scalpel Handle, Round                                 | 00190776028305 |
| 03.520.023     | Caliper, Castroviejo Style, 0mm-40mm                  | 00190776028312 |
| 03.520.024     | Mallet, Small   | 00190776028329 |
| 60.520.001     | Graphic Case/Tissue and Palatal Retractor Instruments | 00190776121204 |

## Tissue and Palatal Retractor Instrumentation System Instructions

### INTENDED USE

- The Tissue and Palatal Retractor Instrumentation System is intended to offer a comprehensive set of surgical instruments for the repair of cleft lip and palate deformities.

### INTENDED USER PROFILE

- Surgical procedures should be performed only by persons having adequate training and familiarity with surgical techniques.
- Consult medical literature relative to techniques, complications and hazards prior to performance of any surgical procedure. Before using the product, all instructions regarding its safety features must be read carefully.

### DEVICE DESCRIPTION

- Surgical instruments comprising fixed assemblies, simple hinged instruments and simple assemblies generally composed of medical grade stainless steels, titanium and aluminum.
- Instrument case and trays may consist of different materials including stainless steels, aluminum and silicone mats.
- Instruments are supplied NON-STERILE and must be inspected, cleaned and sterilized before each use.
- Devices are critical and require terminal sterilization per FDA guidelines and the Spaulding Classification scheme.
- Devices are not implantable.

### WARNINGS

- Avalign recommends thorough manual and automated cleaning of medical devices prior to sterilization. Automated methods alone may not adequately clean devices.
- Devices should be reprocessed as soon as possible following use. Instruments must be cleaned separately from cases and trays.
- All cleaning agent solutions should be replaced frequently before becoming heavily soiled.
- Prior to cleaning, sterilization and use, remove all protective caps carefully. All instruments should be inspected to ensure proper function and condition. Do not use instruments if they do not perform satisfactorily.
- The sterilization methods described have been validated with the devices in predetermined placement locations per the case and tray designs. Areas intended for specific devices shall contain only those devices.
- Proper loading in the desktop gravity displacement unit is critical to the validated sterilization method described.
- Risk of damage – The surgical instrument is a precision device. Careful handling is important for the accurate functioning of the product. Improper external handling can cause product malfunction.
- Use caution when handling sharp instruments to avoid injury.
- Wash the instrument case and trays with an aluminum safe, neutral pH detergent to avoid faded surface colors and deterioration of anodized surfaces.
- If a device is/was used in a patient with, or suspected of having Creutzfeldt-Jakob Disease (CJD), the device cannot be reused and must be destroyed due to the inability to reprocess or sterilize to eliminate the risk of cross-contamination.

### CAUTION

Federal U.S. Law restricts this device to sale, distribution, and use, by, or on order of a physician.

### LIMITATIONS ON REPROCESSING

Repeated processing has minimal effect on these instruments. End of life is normally determined by wear and damage due to use.

### DISCLAIMER

It is the responsibility of the reprocessor to ensure reprocessing is performed using equipment, materials and personnel in the reprocessing facility and achieves the desired result. This requires validation and routine monitoring of the process. Any deviation by the reprocessor from the instructions provided must be properly evaluated for effectiveness and potential adverse consequences.

## Reprocessing Instructions

### TOOLS AND ACCESSORIES

|                 |   |
|-----------------|---|
| Water           | Cold Tap Water (< 20°C / 68°F)<br>Hot Tap Water (> 40°C / 104°F)<br>Deionized (DI) or Reverse Osmosis (RO) Water (ambient)  |
| Cleaning Agents | Neutral Enzymatic Detergent pH 6.0-8.0 i.e. MetriZyme, EndoZime, Enzol  |
| Accessories     | Assorted Sizes of Brushes and/or Pipe Cleaners with Nylon Bristles<br>Sterile Syringes or equivalent<br>Absorbent, Low Lint Disposable Cloths or equivalent<br>Soaking Pans |
| Equipment       | Medical Compressed Air<br>Ultrasonic Cleaner (Sonicator)<br>Automated Washer  |

### POINT-OF-USE AND CONTAINMENT

- 1) Follow health care facility point of use practices. Keep devices moist after use to prevent soil from drying and remove excess soil and debris from all surfaces, crevices, sliding mechanisms, hinged joints, and all other hard-to-clean design features.
- 2) Follow universal precautions and contain devices in closed or covered containers for transport to central supply.

### MANUAL CLEANING

- 3) All devices must be cleaned in the completely open and disassembled (i.e. taken apart) configuration. Disassemble frame blades from palatal retractor frame.
- 4) Prepare neutral pH enzymatic detergent per vendor's directions. Enzol® enzymatic detergent is recommended at a preparation of 1 oz./gallon using lukewarm water.
- 5) Fully immerse device in the prepared detergent per labeling instructions. Allow device to soak for a minimum of 1 minute.
- 6) Actuate all movable parts during the soak time to allow complete penetration of detergent to hard to reach areas.
- 7) Scrub the device, using a soft bristled brush (may also include a syringe and pipe cleaner), paying particular attention to movable parts, crevices, and other hard to reach areas until all visible soil has been removed.
- 8) Prepare neutral pH enzymatic detergent in the sonicator (as per vendor directions) and sonicate the instruments for a minimum of 10 minutes. Note: Enzyme solution shall be changed when it becomes grossly contaminated (bloody and/or turbid).
- 9) Rinse all surfaces and crevices in running reverse osmosis or deionized (RO/DI) water for a minimum of 3 minutes to remove any residual detergent or debris.
- 10) Dry the instrument with a clean, soft cloth. Filtered, compressed air may be used to aid drying.
- 11) Visually examine each instrument for cleanliness. If visible soil remains, repeat cleaning procedure.

### AUTOMATED CLEANING

Note: All devices must be manually pre-cleaned prior to any automated cleaning process, follow steps 1-7. Steps 8-11 are optional but advised.

- 12) Clean devices within a washer/disinfector utilizing the equipment and detergent manufacturers' instructions per the below minimum parameters.

| Phase                | Time (minutes) | Temperature         | Detergent Type & Concentration |
|----------------------|----------------|---------------------|--------------------------------|
| Pre-wash 1           | 02:00          | Cold Tap Water      | N/A                            |
| Enzyme Wash          | 02:00          | Hot Tap Water       | Enzyme Detergent               |
| Rinse 1              | 01:00          | Hot Tap Water       | N/A                            |
| Purified Water Rinse | 00:10          | 146-150°F / 63-66°C | N/A                            |
| Drying               | 15:00          | 194°F / 90°C        | N/A                            |

- 13) Dry excess moisture using an absorbent cloth. Dry any internal areas with filtered, compressed air.
- 14) Visually examine each instrument for cleanliness. If visible soil remains, repeat cleaning procedure.

### DISINFECTION

- Devices must be terminally sterilized (See § Sterilization).
- Aalign instruments are compatible with washer/disinfector time-temperature profiles for thermal disinfection per ISO 15883.
- Load the devices in the washer-disinfector according to the manufacturer's instructions, ensuring that the devices and lumens can drain freely.
- The following automated cycles are examples of validated cycles:

| Phase                | Recirculation Time (min.) | Water Temperature | Water Type  |
|----------------------|---------------------------|-------------------|-------------|
| Thermal Disinfection | 1                         | >90°C (194°F)     | Rl/DO Water |
| Thermal Disinfection | 5                         | >90°C (194°F)     | Rl/DO Water |

### INSPECTION AND FUNCTIONAL TESTING

- Check for smooth movement of hinges. Locking mechanisms should be free of nicks.
- Instruments with broken, cracked, chipped or worn parts should not be used, but should be repaired or replaced immediately.
- Lubricate the instrument before autoclaving with Instra-Lube, or a steam permeable instrument lubricant.

### PACKAGING

- Only FDA cleared sterilization packaging materials should be used by the end user when packaging the devices.
- The end user should consult ANSI/AAMI ST79 or ISO 17665-1 for additional information on steam sterilization.
- **Sterilization Wrap**
  - Prevacuum Cycle: Cases may be wrapped in a standard, medical grade sterilization wrap using a double layer wrap per the AAMI method or equivalent.
  - Gravity Displacement Cycle: Individual trays may be wrapped in a standard, medical grade sterilization wrap using a single layer wrap per the AAMI method or equivalent.
- **Rigid Sterilization Container**
  - For information regarding rigid sterilization containers, please refer to appropriate instructions for use provided by the container manufacturer or contact the manufacturer directly for guidance.

### STERILIZATION

Sterilize with steam. The following are minimum cycles required for steam sterilization of Avalign devices:

#### Double Wrapped Instrument Case:

| Cycle Type | Temperature   | Exposure Time | Pulses | Drying Time |
|------------|---------------|---------------|--------|-------------|
| Prevacuum  | 132°C (270°F) | 4 minutes     | 4      | 25 minutes  |


#### Single Wrapped Instrument Trays:

| Cycle Type           | Temperature   | Cycle Time | Drying Time |
|----------------------|---------------|------------|-------------|
| Gravity Displacement | 134°C (273°F) | 30 minutes | 30 minutes  |


⚠ Gravity Displacement Cycle validated in a Tuttnauer ValueKlave 1730

- Sterilization validation with a 30 minute cycle time resulted in a required minimum exposure time of 14 minutes.
- Metal rack must be inverted 180° to facilitate loading and proper airflow per the efficacy study exposure and dry times. Rack is inverted when the crossbar is at the top of the chamber.
- No more than two trays can be processed at one time in the Tuttnauer unit. Trays must be loaded with tray ledges fully seated on the shelves of the rack. When processing two trays, trays must be loaded on the top and bottom shelves (refer to images below). The tallest tray, Tray 4, can only be loaded on the bottom shelf.
- Trays should be carefully loaded to avoid tearing or damaging the sterilization wrap. For proper airflow and drying, wrapped trays should not touch the chamber wall and should sit above the level of the water reservoir.
- Loading Configuration with inverted rack:

Chamber w/ Unwrapped Trays  
(for visual reference)

Chamber w/ Wrapped Trays

- The operating instructions and guidelines for maximum load configuration of the sterilizer manufacturer should be followed explicitly. The sterilizer must be properly installed, maintained, and calibrated.
- Time and temperature parameters required for sterilization vary according to type of sterilizer, cycle design, and packaging material. It is critical that process parameters be validated for each facility's individual type of sterilization equipment and product load configuration.
- A facility may choose to use different steam sterilization cycles other than the cycle suggested if the facility has properly validated the cycle to ensure adequate steam penetration and contact with the devices for sterilization. Note: rigid sterilization containers cannot be used in gravity steam cycles.
- Water droplets and visible signs of moisture on sterile packaging/wrap or the tape used to secure it may compromise the sterility of the processed loads or be indicative of a sterilization process failure. Visually check outside wrap for dryness. If there are water droplets for visible moisture observed the pack or instrument tray is considered unacceptable. Repackaging and re-sterilize the packages with visible signs of moisture.

#### STORAGE

- After sterilization, instruments should remain in sterilization packaging and be stored in a clean, dry cabinet or storage case.
- Care should be taken when handling devices to avoid damaging the sterile barrier.

#### MAINTENANCE

- Attention: Apply lubricant only on the connecting elements (locking mechanism) and moving parts.
- Discard damaged, worn or non-functional devices.

#### WARRANTY

- All products are guaranteed to be free from defects in material and workmanship at the time of shipping.
- Avalign instruments are reusable and meet AAMI standards for sterilization. All our products are designed and manufactured to meet the highest quality standards. We cannot accept liability for failure of products which have been modified in any way from their original design.

#### NOTICE

- Any serious incident that has occurred in relation to the device(s) should be reported to the Manufacturer, Avalign Technologies Inc, and the competent authority of the EU Member State in which the user and/or patient is established.

#### CONTACT

**Manufactured by:**  
**Avalign Technologies**  
8727 Clinton Park Drive  
Fort Wayne, IN 46825  
1-877-289-1096  
[www.avalign.com](http://www.avalign.com)  
[product.questions@avalign.com](mailto:product.questions@avalign.com)

**Authorized Representative Instruments:**  
**Instrumed GmbH**  
(dba Avalign German Specialty Instruments)  
78 532 Tuttlingen, Unter Buchsteig 3  
Germany Phone: +49 7462 200 49 0

**Distributed in US By:**  
 Synthes USA, LLC  
 1101 Synthes Avenue  
 Monument, CO 80132

**Authorized Representative Case and Tray:**  
 Emergo Europe,  
 Prinsessegracht 20, 2514 AP  
 The Hague, The Netherlands

**Distributed OUS By:**  
 Synthes GMBH  
 Eimattstrasse 3  
 4436 Oberdorf  
 Switzerland

**Symbols Glossary**

| Symbol | Title   | Symbol | Title and Translations  |
|--------|---|--------|---|
|        | Manufacturer & Date of Manufacture                  |        | Caution   |
|        | Lot Number / Batch Code                             |        | Non-Sterile   |
|        | Catalogue Number                                    |        | Federal Law (USA) restricts this device to sale by or on the order of a physician |
|        | Consult Instructions for Use                        |        | Medical Device  |
|        | Authorized Representative in the European Community |        |   |

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