## Disc Preparation Surgical Instruments Instructions for Use

Translations available (click on language to advance to translation):

LANGUAGE
English (EN)
French (Français - FR)
German (Deutsch - DE)
Italian (Italiano - IT)
Portuguese (Português - PT)
Spanish (Español - ES)
Dutch (Nederlands - NL)
Greek (Ελληνικά - EL)
Danish (Dansk – DA)
Finnish (Finnish - FIN)
Swedish (Svenska – SE)
Norwegian (Norsk – NO)
Polish (Polski – PL)
Romanian (Română – RO)
Bulgarian (Български – BG)
Turkish (Türkçe – TR)
Chinese (简体中文 – CN)
Korean (한국어 – KO)

Part Number	Description and Translation	GTIN
48970135	<b>30° Angled Ring Curette</b> Curette arrondie coudée/ Cureta de anillo angulada / Gewinkelte Ringkürette / Curette ad anello angolata/ Ringcurette met hoek / Cureta de anel angulada/ Ξέστρο με δακτυλίδι γωνίας/ Vinklet ringcurette / Kulman rengaskyretti / Chiuretă inel angulată/ Vinklet ringkyrette / Łyzeczka okrągła, kątowa/ Chiuretă inel angulată/ Ёглова пръстеновидна кюрета / Açılı Halka Küret / 成角环式刮匙 / 앵글드 링큐렛	00190776138233
48970015	<b>Straight 4mm Cup Curette</b> Curette coupelle droite/ Cureta de copa recta/ Gerade Kürette mit Kappe/ Curette a cucchiaio retta/ Rechte bekerurette/ Cureta de Volkmann/ Ξέστρο με ευθεία κόψα / Lige kopcurette/ Suora kuppikyretti/ Chiuretă cupă dreaptă / Rett kopkyrette/ Łyzeczka z główką, prosta/ Chiuretă cupă dreaptă/ Права кюрета с чашка/ Düz Oval Uçlu Küret / 直型 杯式刮匙 / 일자형 컵 큐렛	00190776140557
48970017	<b>45° 4mm Cup Curette</b> Curette coupelle / Cureta de copa angulada/ Gewinkelte Kürette mit Kappe/ Curette a cucchiaio/ Bekercurette/ Cureta de Volkmann/ Ξέστρο με ευθεία κόψα / Kopcurette/ Kuppikyretti/ Chiuretă cupă/ Kopkyrette/ Łyzeczka z główką/ Chiuretă cupă/ кюрета с чашка/ Oval Uçlu Küret/ 杯式刮匙 / 컵 큐렛	00190776140571
48970016	<b>Straight 7mm Cup Curette</b> Curette coupelle droite/ Cureta de copa recta/ Gerade Kürette mit Kappe/ Curette a cucchiaio retta/ Rechte bekerurette/ Cureta de Volkmann/ Ξέστρο με ευθεία κόψα / Lige kopcurette/ Suora kuppikyretti/ Chiuretă cupă dreaptă / Rett kopkyrette/ Łyzeczka z główką, prosta/ Chiuretă cupă dreaptă/ Права кюрета с чашка/ Düz Oval Uçlu Küret / 直型 杯式刮匙 / 일자형 컵 큐렛	00190776140564
48970018	<b>45° 7mm Cup Curette</b> Curette coupelle / Cureta de copa angulada/ Gewinkelte Kürette mit Kappe/ Curette a cucchiaio/ Bekercurette/ Cureta de Volkmann/ Ξέστρο με ευθεία κόψα / Kopcurette/ Kuppikyretti/ Chiuretă cupă/ Kopkyrette/ Łyzeczka z główką/ Chiuretă cupă/ кюрета с чашка/ Oval Uçlu Küret/ 杯式刮匙 / 컵 큐렛	00190776140588
48970021	<b>Straight Osteotome, 10mm</b> Ostéotome droit/ Osteótomo recto/ Gerades Osteotom/ Osteotomo retto/ Rechte osteotom/ Osteótomo recto/ Ευθύ εργαλείο οστεοτομής/ Lige osteotom/ Suora osteotomi/ Osteotom drept/ Rett osteotom/ Osteotom prosty/ Osteotom drept/ Прав остеотом/ Düz Osteotom/ 直型骨齒 / 일자형 뼈절단기	00190776140618
48970022	<b>Curved Osteotome, 10mm</b> Ostéotome incurvé/ Osteótomo curvado/ Gekrümmtes Osteotom/ Osteotomo curvo/ Gebogen osteotom/ Osteótomo curvo/ Κυρτό εργαλείο οστεοτομής/ Buget osteotom/ Kaareva osteotomi/ Osteotom curbat/ Kurvet osteotom/ Osteotom zagięty/ Osteotom curbat/ Извит остеотом/ Eğri Osteotom/ 弯型骨齒 / 곡선형 뼈절단기	00190776140625
48970129	<b>Pituitary 4mm, Non-Bayonnetted</b> Pituitaire, non à baïonnette/ Piyuitario, sin forma de bayoneta/ Hypophysen-Rongeur, ohne Bajonett/ Pituitaria, senza baionetta/ Hypofyse, zonder bayonet/ Pituitária sem baioneta/ Υπόφυση, χωρίς μπαγιονέτ/ Hypofyse, uden bajonet/ Pituitaarinen, pistimetön/ Pituitară, fără baionetă/ Pituitary, ikke-bajonett/ Odgryzacz typu Pituitary, niebagnetowy/ Pituitară, fără baionetă/ За хипофиза, без шип/ Pituiter, Bayonetsiz/ 垂体, 非刺刀式/ 뇌하수체용, 비총검형	00190776141028
48970130	<b>Pituitary 6mm, Non-Bayonnetted</b> Pituitaire, non à baïonnette/ Piyuitario, sin forma de bayoneta/ Hypophysen-Rongeur, ohne Bajonett/ Pituitaria, senza baionetta/ Hypofyse, zonder bayonet/ Pituitária sem baioneta/ Υπόφυση, χωρίς μπαγιονέτ/ Hypofyse, uden bajonet/ Pituitaarinen, pistimetön/ Pituitară, fără baionetă/ Pituitary, ikke-bajonett/ Odgryzacz typu Pituitary, niebagnetowy/ Pituitară, fără baionetă/ За хипофиза, без шип/ Pituiter, Bayonetsiz/ 垂体, 非刺刀式/ 뇌하수체용, 비총검형	00190776141035
48970131	<b>Pituitary Up-biting 4mm, Non-Bayonnetted</b> Pituitaire à mors supérieurs, non à baïonnette/ Piyuitario, sin forma de bayoneta/ Hypophysen-Rongeur, oben, ohne Bajonett/ Pituitaria Upbiting, senza baionetta/ Hypofyse omhoog happened, zonder bayonet/ Pituitária de preensão superior sem baioneta/ Υπόφυση άνω σκαλίσματος, χωρίς μπαγιονέτ/ Hypofyse m. opadvendt bid, uden bajonet/ Pituitaarinen, ylös pureva, pistimetön/ Pituitară cu mușcătura sus, fără baionetă/ Pituitary opp-bitende, ikke-bajonett/ Odgryzacz typu Pituitary odcinanie górne, niebagnetowy/ Pituitară cu mușcătura sus, fără baionetă/ За хипофиза, с отвор нагоре, без шип/ Pituiter Üst Isırma, Bayonetsiz/ 垂体上咬式, 非刺刀式/ 뇌하수체용 Up-biting, 비총검형	00190776141042

Part Number	Description and Translation	GTIN
48970132	<b>Kerrison 3mm, Non-Bayoneted</b> Kerrison, non à baïonnette/ Piyuitario, mordida hacia arriba, sin forma de bayoneta/ Kerrison-Rongeur, ohne Bajonett/ Kerrison, senza baionetta/ Kerrison, zonder bayonet/ Kerrison sem baioneta/ Kerrison, χωρίς μπαγιονέτ/ Kerrison, uden bajonet/ Kerrison, pistimetön/ Kerrison, fără baionetă/ Kerrison, ikke-bajonett/ Kerrison, niebagnetowy/ Kerrison, fără baionetă/ Kerrison, без шип/ Kerrison, Bayonetsiz/ Kerrison, 非刺刀式/ Kerrison, 비총검형	00190776141059
48970133	<b>Kerrison 5mm, Non-Bayoneted</b> Kerrison, non à baïonnette/ Piyuitario, mordida hacia arriba, sin forma de bayoneta/ Kerrison-Rongeur, ohne Bajonett/ Kerrison, senza baionetta/ Kerrison, zonder bayonet/ Kerrison sem baioneta/ Kerrison, χωρίς μπαγιονέτ/ Kerrison, uden bajonet/ Kerrison, pistimetön/ Kerrison, fără baionetă/ Kerrison, ikke-bajonett/ Kerrison, niebagnetowy/ Kerrison, fără baionetă/ Kerrison, без шип/ Kerrison, Bayonetsiz/ Kerrison, 非刺刀式/ Kerrison, 비총검형	00190776141066
48970056	<b>Single Sided Rasp</b> Râpe simple face/ Lima de un lado/ Einseitige Raspel/ Raspa monolaterale/ Eenzijdige rasp/ Raspador de face simples/ Ράσπα μονής πλευράς/ Enkeltsidet rasp/ Yksipuolinen raspi/ Raşpel cu o faţă/ Enkeltsidet rasp/ Skrobaczka jednostronna/ Raşpel cu o faţă/ Одностронна пила/ Tek Taraflı Raspa/ 单侧骨锉/ 단일 면 뼈막막리기	00190776140960
48970057	<b>Double Sided Rasp</b> Râpe double face/ Lima de dos lados/ Doppelseitige Raspel/ Raspa bilaterale/ Dubbelzijdige rasp/ Raspador de dupla face/ Ράσπα διπλής πλευράς/ Dobbelsidet rasp/ Kaksipuolinen raspi/ Raşpel cu două feţe/ Dobbelsidet rasp/ Skrobaczka dwustronna/ Raşpel cu două feţe/ Двустронна пила/ Çift Taraflı Raspa/ 双侧骨锉/ 이중 면 뼈막막리기	00190776140977
48970058	<b>Penfield 4mm, Bayoneted</b> Penfield, à baïonnette / Penfield, en bayoneta / Penfield-Dissektor, mit Bajonett / Penfield, a baionetta / Penfield, met bayonet / Penfield com baioneta/ Penfield, με μπαγιονέτ / Penfield, med bajonet / Penfield, pistimellinen / Penfield, tip baionetă / Penfield, bajonett / Penfield, bagnetowy / Penfield, tip baionetă / Penfield, с шип / Penfield, Bayonetli / Penfield, 刺刀式 / Penfield, 총검형	00190776140984
48970059	<b>Penfield Reverse 4mm, Bayoneted</b> Penfield inversé, à baïonnette / Con cambio Penfield, en bayoneta / Penfield-Rückwärts-Dissektor, mit Bajonett / Penfield inverso, a baionetta / Penfield omgekeerd, met bayonet / Penfield Inverso com baioneta/ Ανάστροφο Penfield, με μπαγιονέτ / Penfield, omvendt, med bajonet / Penfield, käänteinen, pistimellinen / Penfield inversat, tip baionetă / Penfield revers, bajonett / Penfield odwrócony, bagnetowy / Penfield inversat, tip baionetă / Penfield обратен, с шип / Penfield Ters, Bayonetli / Penfield 翻转 · 刺刀式 / Penfield 후방용, 총검형	00190776140991
48970062	<b>Knife Handle, Bayoneted</b> Manche de couteau, à baïonnette / Mango de cuchilla, en bayoneta / Messergriff, mit Bajonett / Impugnatura coltello, a baionetta / Messenheft, met bayonet / Cabo de bisturi com baioneta / Λαβή μαχαριού, με μπαγιονέτ / Knivhåndtag, med bajonet / Veitsenkahva, pistimellinen / Mâner cuțit, tip baionetă / Knivhåndtak, bajonett / Trzonek do ostrza, bagnetowy / Mâner cuțit, tip baionetă / Дръжка за нож, с шип / Viçak Sarı, Bayonetli / 刀柄, 刺刀式 / 나이프 손잡이, 총검형	00190776141004
48019301	<b>Slap Hammer</b> Impacteur marteau/ Martillo de dos caras/ Ausschläger/ Martello percussore/ Slaghamer/ Martelo cirúrgico/ Σφυρί / Kirurgisk hammer/ Slap-vasara/ Ciocan chirurgical/ Kirurgisk hammer/ Wybijak/ Ciocan chirurgical/ Чук/ Cerrahi Çekiç/ 手术锤/ 슬랩 해머	00190776138271
48970098	<b>Mallet</b> Maillet / Mazo / Schlägel / Mazzetta / Hamer / Malhete / Ξύλινο σφυρί / Hammer / Vasara / Ciocan / Klubbe / Pobijak / Ciocan / Чукче / Çekiç / 槌 / 망치	00190776141318
48970067	<b>Frazier Suction 11Fr</b> Canule d'aspiration type Frazier/ Tubo de aspiración de Frazier/ Frazier-Ansaugschlauch/ Aspirazione di Frazier/ Frazier afzuiging/ Tubo de aspiração Frazier/ Βεντούζα Frazier/ Frazier sugearordning/ Frazier-imu/ Tub de aspirație Frazier/ Frazier sug/ Rurka odsysająca Frazier/ Tub de aspirație Frazier/ Аспирационна тръба Frazier/ Frazier Emme/ Frazier 吸引 / 프레이저 흡입기	00190776141011
48970127	<b>Frazier Suction 14Fr</b> Canule d'aspiration type Frazier/ Tubo de aspiración de Frazier/ Frazier-Ansaugschlauch/ Aspirazione di Frazier/ Frazier afzuiging/ Tubo de aspiração Frazier/ Βεντούζα Frazier/ Frazier sugearordning/ Frazier-imu/ Tub de aspirație Frazier/ Frazier sug/ Rurka odsysająca Frazier/ Tub de aspirație Frazier/ Аспирационна тръба Frazier/ Frazier Emme/ Frazier 吸引 / 프레이저 흡입기	00190776141561

## Label Glossary

Symbol	Title and Translations
	Manufacturer/ Fabricant/ Hersteller/ Fabbricante/ Fabricante/ Fabricante/ Fabrikant/ Κατασκευαστής/ Producent/ Valmistaja/ Tillverkare/ Produzent/ Producent/ Fabricant/ Производител/ Üretici/ 制造商/ 제조사
	Lot Number / Batch Code/ Numéro de lot/Code de lot/ Chargennummer/ Numero di lotto/codice partita/ Número de lote / Código de série/ Número de lote / Batch Code/ Lotnummer / Partijcode/ Αριθμός Παρτίδας / Κωδικός Παρτίδας/ Serienummer / partikode/ Eränumero/eräkoodi/ Satsnummer/batchkod/ Lotnummer / Batchkode/ Nr partii / Kod partii/ Număr lot/Cod serie/ Партида № / Код на партида/ Lot Numarası / Parti Kodu/ 批号/批次代码/ 로트 번호/배치 코드
	Catalogue Number/ Référence/ Katalognummer/ Numero catalogo/ Número de referência/ Número de catálogo/ Catalogusnummer/ Αριθμός Καταλόγου/ Katalognummer/ Luettelonumero/ Katalognummer/ Katalognummer/ Numer katalogowy/ Număr de catalog/ Καταλογен номер/ Katalog Numarası/ 目录号/ 카탈로그 번호
	Consult Instructions for Use/ Consulter le mode d'emploi/ Gebrauchsanweisung beachten/ Consultare le istruzioni per l'uso/ Consultar as Instruções de Utilização/ Consultar las instrucciones de uso/ Raadpleeg de gebruiksaanwijzingen/ Συμβουλευθείτε τις Οδηγίες Χρήσης/ Læs brugsanvisningen/ Tutustu käyttöohjeisiin/ Bruksanvisning/ Se i bruksanvisningen/ Sprawdzić w instrukcji użycia/ Consultați instrucțiunile de utilizare/ Консултирайте се с указанията за употреба/ Kullanım Talimatlarına Başvurun/ 参阅使用说明/ 사용 설명서 참조
	Medical Device/ Dispositif médical/ Medizinprodukt/ Dispositivo medico/ Dispositivo Médico/ Dispositivo médico/ Medisch instrument/ Ιατρική συσκευή/ Medicinsk udstyr/ Lääkinnällinen laite/ Medicinsk enhet/ Medisinsk enhet/ Wyrób medyczny/ Dispozitiv medical/ Медицинско изделие/ Tibbi Cihaz/ 医疗装置/ 의료기기
	Caution/ Attention/ Achtung/ Attenzione/ Aviso/ Precaución/ Waarschuwing/ Προσοχή/ Forsigtig/ Varoitus/ Försiktighet/ Forsiktig/ Uwaga/ Atenție/ Внимание/ Dikkat/ 注意/ 주의사항(Caution)
	Non-Sterile/ Non stérile/ Unsteril/ Non sterile/ Não estéril/ No estéril/ Niet steriel/ Μη αποστειρωμένο/ Ikke-steril/ Ei-steriili/ Icke-steril/ Ikke-steril/ Niesterylny/ Nesteril/ Не стерилно/ Steril Değildir/ 非无菌/ 비멸균(Non-Sterile)
	Federal Law (USA) restricts this device to sale by or on the order of a physician/ La législation fédérale américaine limite la vente de ce dispositif à un médecin ou sur ordre d'un médecin/ Nach US-amerikanischem Bundesgesetz darf dieses Gerät nur von einem zugelassenen Arzt oder auf ärztliche Anordnung verkauft werden/ La legge federale (USA) limita la vendita di questo dispositivo da parte o su prescrizione di un medico/ A Lei Federal (EUA) restringe a venda deste dispositivo a médicos credenciados ou mediante receita médica/ La Ley Federal de Estados Unidos limita la venta de este dispositivo a un médico o por prescripción de éste/ Volgens de federale wetgeving (VS) mag dit instrument uitsluitend door of in opdracht van een arts worden gekocht/ Ο Ομοσπονδιακός Νόμος (Η.Π.Α.) επιτρέπει την πώληση αυτής της συσκευής μόνο σε ή κατόπιν εντολής ιατρού/ Ifølge amerikansk lov må dette udstyr kun sælges af eller efter ordination fra en læge/ Yhdysvaltain liittovaltion laki rajoittaa tämän laitteen myynnin tapahtuvan ainoastaan lääkärille tai lääkärin määräyksestä/ Enligt federal lag (USA) får denna enhet endast säljas av legitimerad läkare eller efter läkares föreskrift./ Føderal lov (USA) begrænser denne innretningen til salg av eller etter ordre fra lege./ Zgodnie z prawem federalnym (USA) wyrób ten może być sprzedawany lub zamawiany wyłącznie przez lekarzy/ Legea federală (SUA) restricționează vânzarea acestui dispozitiv numai de către un medic sau în baza prescripției acestuia/ Федералното законодателство (на САЩ) ограничава продажбата на това изделие от или по поръчка на лекар/ ABD federal kanunlarına göre bu cihaz sadece bir doktor tarafından ya da doktorun siparişiyle satılabilir/ 美国联邦法律规定，本装置只能由有执照的医生销售或凭医嘱销售。/ 연방법(USA)상 이 기기는 의사에게만 판매 또는 주문되도록 제한되어 있습니다.
	Authorized Representative in the European Community/ Représentant agréé dans la Communauté européenne/ Autorisierte Vertretung in der Europäischen Gemeinschaft/ Rappresentante autorizzato nella Comunità europea/ Representante Autorizado na União Europeia/ Representante autorizado en la Comunidad Europea/ Geautoriseerd vertegenwoordiger in de Europese Gemeenschap/ Εξουσιοδοτημένος Αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα/ Autoriseret repræsentant i Den Europæiske Union/ Valtuutettu edustaja Euroopan yhteisössä/ Auktoriserad EU-representant/ Autorisert representant i Det europeiske fellesskap/ Autoryzowany przedstawiciel we Wspólnocie Europejskiej/ Reprezentant autorizat pentru Comunitatea Europeană/ Упълномощен представител в Европейската общност/ Avrupa Birliği Yetkili Temsilcisi/ 欧共同体授权代表/ 유럽 연합 내의 위임 대표

## Disc Preparation Surgical Instruments Instructions for Use

### INTENDED USE

- The Disc Preparation Surgical Instruments are intended to offer a comprehensive set of surgical instruments to prepare the intervertebral disc space for interbody spinal fusion.

### INTENDED USER PROFILE

- Surgical procedures should be performed only by persons having adequate training and familiarity with surgical techniques.
- Consult medical literature relative to techniques, complications and hazards prior to performance of any surgical procedure. Before using the product, all instructions regarding its safety features must be read carefully.

### DEVICE DESCRIPTION

- Surgical instruments comprising fixed assemblies, simple hinged instruments and simple assemblies generally composed of medical grade stainless steels, titanium, aluminum and silicone rubber.
- Devices are supplied NON-STERILE and must be inspected, cleaned and sterilized before each use.
- Devices are critical and require terminal sterilization.
- Devices are not implantable.
- Frazier Suction Tubes interface with standard surgical tubing sizes over the barbed fitting: 6mm, ¼" (6.35mm), 7mm, 9/32" (7.14mm), 5/16" (7.94mm), and 8mm.

### WARNINGS

- Avalign recommends thorough manual and automated cleaning of medical devices prior to sterilization. Automated methods alone may not adequately clean devices.
- Devices should be reprocessed as soon as possible following use. Instruments must be cleaned separately from cases and trays.
- **For the Frazier Suction Tubes, disassemble and remove stylet before cleaning.**
- All cleaning agent solutions should be replaced frequently before becoming heavily soiled.
- Prior to cleaning, sterilization and use, remove all protective caps carefully. All instruments should be inspected to ensure proper function and condition. Do not use instruments if they do not perform satisfactorily.
- **The sterilization methods described have been validated with the devices in predetermined placement locations per the case and tray designs. Areas intended for specific devices shall contain only those devices.**
- Proper installation and maintenance instructions for interfacing suction equipment should be followed to achieve desired vacuum pressure and flow rate.
- If vacuum pressure or flow rate are diminished, verify tubing is installed over at least one of the suction instrument's barbed connection features. Utilize stylet to remove any material lodged in the suction cannula.
- Risk of damage – The surgical instruments are precision devices. Careful handling is important for the accurate functioning of the devices. Improper external handling can cause the devices to malfunction.
- Use caution when handling sharp instruments to avoid injury.
- If a device is/was used in a patient with, or suspected of having Creutzfeldt-Jakob Disease (CJD), the device cannot be reused and must be destroyed due to the inability to reprocess or sterilize to eliminate the risk of cross-contamination.

### CAUTION

Federal U.S. Law restricts this device to sale, distribution, and use, by, or on order of a physician.

### LIMITATIONS ON REPROCESSING

Repeated processing has minimal effect on these instruments. End of life is normally determined by wear and damage due to use.

### DISCLAIMER

It is the responsibility of the reprocessor to ensure reprocessing is performed using equipment, materials and personnel in the reprocessing facility and achieves the desired result. This requires validation and routine monitoring of the process. Any

deviation by the reprocessor from the instructions provided must be properly evaluated for effectiveness and potential adverse consequences.

## Reprocessing Instructions

### TOOLS AND ACCESSORIES

Water	Cold Tap Water (< 20°C / 68°F) Warm Water (38°- 49°C / 100°- 120°F) Hot Tap Water (> 40°C / 104°F) Deionized (DI) or Reverse Osmosis (RO) Water (ambient)
Cleaning Agents	Neutral Enzymatic Detergent pH 6.0-8.0 i.e. MetriZyme, EndoZime, Enzol
Accessories	Assorted Sizes of Brushes and/or Pipe Cleaners with Nylon Bristles Sterile Syringes or equivalent Absorbent, Low Lint Disposable Cloths or equivalent Soaking Pans
Equipment	Medical Compressed Air Ultrasonic Cleaner Automated Washer

### POINT-OF-USE AND CONTAINMENT

- 1) Follow health care facility point of use practices. Keep devices moist after use to prevent soil from drying and remove excess soil and debris from all lumens, surfaces, crevices, sliding mechanisms, hinged joints, and all other hard-to-clean design features.
- 2) Suction or flush lumens with a cleaning solution immediately after use.
- 3) Follow universal precautions and contain devices in closed or covered containers for transport to central supply.

### MANUAL CLEANING

- 4) Disassemble all devices as warranted per manufacturer's instructions.
- 5) Rinse devices under cold running tap water for a minimum of 3 minutes while wiping off residual soil or debris. Actuate moveable mechanisms and flush all lumens, cracks and/or crevices while rinsing.
- 6) Prepare an enzymatic cleaning solution per manufacturer's instructions including dilution/concentration, water quality and temperature. Immerse devices and soak for a minimum of 10 minutes. While in the solution, use a soft, bristle brush to remove all traces of blood and debris from the device, paying close attention to threads, crevices, seams, and any hard to reach areas.
  - a. If the device has sliding mechanisms or hinged joints, actuate the device while scrubbing to remove trapped soil.
  - b. If the device contains a lumen, use a tight-fitting nylon brush or pipe cleaner while pushing in and out with a twisting motion to facilitate removal of debris; ensure the full diameter and depth of the lumen is accessed. Flush the lumen, three times minimum, with a syringe containing a minimum solution of 60mL.
- 7) Remove devices and rinse/agitate in cold tap water for a minimum of 3 minutes. Actuate moveable mechanisms and flush all lumens, cracks and/or crevices while rinsing.
- 8) Prepare a neutral detergent cleaning solution per manufacturer's instructions including dilution/concentration, water quality and temperature. Immerse devices and soak for a minimum of 5 minutes. While in the solution, use a soft, bristle brush to remove all traces of blood and debris from the device, paying close attention to threads, crevices, seams, and any hard to reach areas.
  - a. If the device has sliding mechanisms or hinged joints, actuate the device while scrubbing to remove trapped soil.
  - b. If the device contains a lumen, use a tight-fitting nylon brush or pipe cleaner while pushing in and out with a twisting motion to facilitate removal of debris; ensure the full diameter and depth of the lumen is accessed. Flush the lumen, three times minimum, with a syringe containing a minimum solution of 60mL.
- 9) Remove devices and rinse/agitate in cold tap water for a minimum of 3 minutes. Actuate moveable mechanisms and flush all lumens, cracks and/or crevices while rinsing.
- 10) Prepare an enzymatic cleaning solution using hot water per manufacturer's recommendations in an ultrasonic unit. Sonicate the devices for a minimum of 15 minutes using a minimum frequency of 40 kHz. It is recommended to use an ultrasonic unit with flushing attachments. Devices with lumens should be flushed with cleaning solution under the surface of the solution to ensure adequate perfusion of channels.
- 11) Remove devices and rinse/agitate in ambient DI/RO water for a minimum of 4 minutes. Actuate moveable mechanisms and flush all lumens, cracks and/or crevices while rinsing. Flush internal lumens a minimum of 3 times with RO/DI water (minimum of 15mL) using an appropriately sized syringe. If available, use flush ports for flushing.
- 12) Dry the device using an absorbent cloth. Dry any internal areas with filtered, compressed air.

- 13) Visually inspect the device for soil under magnification including all actuating mechanisms, cracks, crevices, and lumens. If not visibly clean, repeat steps 4-13.
- 14) Submerge device in 2-3% hydrogen peroxide. The appearance of bubbles confirms the presence of hemoglobin. Repeat steps 5-14 if bubbles appear. Adequately rinse device with DI/RO water.

#### AUTOMATED CLEANING

Note: All devices must be manually pre-cleaned prior to any automated cleaning process, follow steps 1-9. Steps 10-14 are optional but advised.

- 15) Transfer the devices to an automatic washer/disinfector for processing per the below minimum parameters.

Phase	Time (minutes)	Temperature	Detergent Type & Concentration
Pre-wash 1	02:00	Cold Tap Water	N/A
Enzyme Wash	02:00	Hot Tap Water	Enzyme Detergent
Wash 1	02:00	63°C / 146°F	Neutral Detergent
Rinse 1	02:00	Hot Tap Water	N/A
Purified Water Rinse	02:00	146°F / 63°C	N/A
Drying	15:00	194°F / 90°C	N/A

- 16) Dry excess moisture using an absorbent cloth. Dry any internal areas with filtered, compressed air.
- 17) Visually inspect the device for soil under magnification including all actuating mechanisms, cracks, crevices and lumens. If not visibly clean, repeat steps 4-9, 15-17.
- 18) Submerge device in 2-3% hydrogen peroxide. The appearance of bubbles confirms the presence of hemoglobin. Repeat steps 5-9, 15-18 if bubbles appear. Adequately rinse device with DI/RO water.

#### DISINFECTION

- Devices must be terminally sterilized (See § Sterilization).
- Aalign instruments are compatible with washer/disinfector time-temperature profiles for thermal disinfection per ISO 15883.

#### INSPECTION AND FUNCTIONAL TESTING

- Visually inspect devices for damage or wear, including sharp edges. Instruments with broken, cracked, chipped or worn features, should not be used, but should be replaced immediately.
- Verify device interfaces (junctions and threads) continue to function as intended without complications.
- Check for smooth movement of hinges. Locking mechanisms should be free of nicks.
- Lubricate before autoclaving with Instra-Lube, or a steam permeable instrument lubricant.

#### PACKAGING

- Only FDA cleared sterilization packaging materials should be used by the end user when packaging the devices.
- The end user should consult ANSI/AAMI ST79 or ISO 17665-1 for additional information on steam sterilization.
- **Sterilization Wrap**
  - Cases may be wrapped in a standard, medical grade sterilization wrap using the AAMI double wrap method or equivalent.
- **Rigid Sterilization Container**
  - For information regarding rigid sterilization containers, please refer to appropriate instructions for use provided by the container manufacturer or contact the manufacturer directly for guidance.

#### STERILIZATION

Sterilize with steam. The following are minimum cycles required for steam sterilization of Aalign devices:

##### 1. US Sterilization Wraps:

Cycle Type	Temperature	Exposure Time	Pulses	Drying Time
Prevacuum	132°C (270°F)	4 minutes	3	30 minutes
Gravity Displacement	132°C (270°F)	15 minutes	-	30 minutes

##### 2. OUS Sterilization Wraps:

Cycle Type	Temperature	Exposure Time	Pulses	Drying Time
Prevacuum	134°C (273°F)	3 minutes	3	30 minutes

##### 3. US Rigid Containers:

Cycle Type	Temperature	Exposure Time	Pulses	Drying Time
Prevacuum	132°C (270°F)	4 minutes	3	30 minutes

#### 4. OUS Rigid Containers:

Cycle Type	Temperature	Exposure Time	Pulses	Drying Time
Prevacuum	134°C (273°F)	3 minutes	3	30 minutes

- The operating instructions and guidelines for maximum load configuration of the sterilizer manufacturer should be followed explicitly. The sterilizer must be properly installed, maintained, and calibrated.
- Time and temperature parameters required for sterilization vary according to type of sterilizer, cycle design, and packaging material. It is critical that process parameters be validated for each facility's individual type of sterilization equipment and product load configuration.
- A facility may choose to use different steam sterilization cycles other than the cycle suggested if the facility has properly validated the cycle to ensure adequate steam penetration and contact with the devices for sterilization. Note: rigid sterilization containers cannot be used in gravity steam cycles.
- Water droplets and visible signs of moisture on sterile packaging or the tape used to secure it, may compromise sterility of processed loads or be indicative of a sterilization process failure. Visually check outside wrapper for dryness. If there are water droplets or visible moisture on the exterior of the package or on the tape used to secure it, the pack or instrument tray is considered unacceptable. Repackage and re-sterilize sterilization packages with visible signs of moisture.

#### STORAGE

- After sterilization, instruments should remain in sterilization packaging and be stored in a clean, dry cabinet or storage case.
- Care should be taken when handling devices to avoid damaging the sterile barrier.

#### MAINTENANCE

- Attention: Apply lubricant only on the connecting elements (locking mechanism) and moving parts.
- Discard damaged, worn or non-functional devices.

#### WARRANTY

- All products are guaranteed to be free from defects in material and workmanship at the time of shipping.
- Aalign instruments are reusable and meet AAMI standards for sterilization. All our products are designed and manufactured to meet the highest quality standards. We cannot accept liability for failure of products which have been modified in any way from their original design.

#### CONTACT

**Manufactured by:**

**Aalign Technologies**

8727 Clinton Park Drive

Fort Wayne, IN 46825

1-877-289-1096

[www.avalign.com](http://www.avalign.com)

[product.questions@avalign.com](mailto:product.questions@avalign.com)

**Authorized Representative:**

**Instrumed GmbH**

(dba Aalign German Specialty Instruments)

78 532 Tuttlingen, Unter Buchsteig 3

Germany Phone: +49 7462 200 49 0

CE mark is valid if printed on product label.

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