Depuy Synthes ATP/Lateral Discectomy Surgical Instrumentation System Instructions

Translations available (click on language to advance to translation):

LANGUAGE
English (EN)
French (Français - FR)
German (Deutsch - DE)
Italian (Italiano - IT)
Portuguese (Português - PT)
Spanish (Español - ES)
Bulgarian (Български - BG)
Chinese (简体中文 - CN)
Croatian (Hrvatski - HR)
Czech (Česky - CS)
Danish (Dansk - DA)
Dutch (Nederlands - NL)
Estonian (Eesti - ET)
Finnish (Finnish - FIN)
Greek (Ελληνικά - EL)
Hungarian (Magyar - HU)
Indonesian (Indonesian - ID)
Japanese (日本語 - JP)
Kazakh (Қазақша - KK)
Korean (한국어 - KO)
Latvian (Latviski - LV)
Lithuanian (Lietuvių k. - LT)
Norwegian (Norsk - NO)
Polish (Polski - PL)
Romanian (Română - RO)
Russian (Русский - RU)
Serbian (Srpski - SR)
Slovak (Slovenčina - SK)
Slovenian (Slovenščina - SL)
Swedish (Svenska - SE)
Thai (ภาษาไทย - TH)
Turkish (Türkçe - TR)
Vietnamese (Tiếng Việt - VI)

Part List

Catalog Number	Description	UDI
287301001	ATP LAT DISCECTOMY LID	00190776152017
287301010	ATP LAT DISCECTOMY CASE 1	00190776151973
287301020	ATP LAT DISCECTOMY CASE 2	00190776151980
287301030	ATP ANGLED INSTRUMENTS CASE	00190776151997
287301040	ATP LAT DISCECTOMY AUX CASE	00190776152000
289215233	TRAY FOR RONGEURS	00190776155384
289215236	TRAY F/STRAIGHT/SGL BEND INSTR	00190776155391
289215241	TRAY F/DINGO/DBL BEND INSTR	00190776155407
287303116	LAT SPREADER 16x5MM	00190776152123
287303206	LAT 6MM ANGLED SPREADER	00190776152918
287304400	LAT 4MM ATP RASP	00190776152369
287304401	LATERAL 4MM STRAIGHT RASP	00190776152208
287304402	LATERAL 4MM ANGLED RASP	00190776152352
287352240	LAT 2MM 40 DEG KERRISON RONGUER	00190776152031
287352440	LAT 4MM 40 DEG KERRISON RONGUER	00190776152048
287353201	LAT 2MM STRAIGHT PITUITARY	00190776152079
287353202	LAT 2MM ANGLED PITUITARY RIGHT	00190776152390
287353203	LAT 2MM ANGLED PITUITARY LEFT	00190776152406
287353204	LAT 2MM UP-ANGLED PITUITARY	00190776152109
287353206	LAT 2MM ANGLED UP-ANGLED PITUITARY RIGHT	00190776152703
287353207	LAT 2MM ANGLED UP-ANGLED PITUITARY LEFT	00190776152710
287353212	LAT 2MM ANGLED PITUITARY RIGHT	00190776151140
287353213	LAT 2MM ANGLED PITUITARY LEFT	00190776151164
287353216	LAT 2MM ANG UP-ANG PITUITARY R	00190776151195
287353217	LAT 2MM ANG UP-ANG PITUITARY L	00190776151225
287353401	LAT 4MM STRAIGHT PITUITARY	00190776152086
287353402	LAT 4MM ANGLED PITUITARY RIGHT	00190776152291
287353403	LAT 4MM ANGLED PITUITARY LEFT	00190776152307
287353404	LAT 4MM UP-ANGLED PITUITARY	00190776152116
287353406	LAT 4MM ANGLED UP-ANGLED PITUITARY RIGHT	00190776152413
287353407	LAT 4MM ANGLED UP-ANGLED PITUITARY LEFT	00190776152420
287353416	LAT 4MM ANG UP-ANG PITUITARY R	00190776151188
287353417	LAT 4MM ANG UP-ANG PITUITARY L	00190776151218
287353601	LAT 6MM STRAIGHT PITUITARY	00190776152093
287353602	LAT 6MM ANGLED PITUITARY, RIGHT	00190776152673
287353603	LAT 6MM ANGLED PITUITARY, LEFT	00190776152680
287353604	LAT 6MM UP-ANGLED PITUITARY	00190776152697
287353606	LAT 6MM ANGLED UP-ANGLED PITUITARY, RIGHT	00190776152727
287353607	LAT 6MM ANGLED UP-ANGLED PITUITARY, LEFT	00190776152734
287353612	LAT 6MM ANGLED PITUITARY RIGHT	00190776151157
287353613	LAT 6MM ANGLED PITUITARY, LEFT	00190776151171
287353614	LAT 6MM UP-ANGLED PITUITARY	00190776151249
287353616	LAT 6MM ANG UP-ANG PITUITARY R	00190776151201
287353617	LAT 6MM ANG UP-ANG PITUITARY L	00190776151232
287354401	LAT 4MM STRAIGHT CURETTE	00190776152161
287354411	LAT 4MM ANGLED CURETTE, RIGHT	00190776152314
287354412	LAT 4MM ANGLED CURETTE, LEFT	00190776152321
287354421	LAT 4MM ATP CUP CURETTE, CRAN	00190776152437
287354422	LAT 4MM ATP CUP CURETTE, CAUD	00190776152444
287354431	LAT 4MM UP-BITING CURETTE	00190776152185
287354441	LAT 4MM ANGLED UP-BITING CURETTE, RIGHT	00190776152451
287354442	LAT 4MM ANGLED UP-BITING CURETTE, LEFT	00190776152468
287354451	LAT 4MM UP-BITING ATP CURETTE	00190776152819
287354453	LAT 4MM DOWN-BITING CURETTE	00190776152475
287354454	LAT 4MM ANGLED DOWN-BITING CURETTE, RIGHT	00190776152833
287354455	LAT 4MM ANGLED DOWN-BITING CURETTE, LEFT	00190776152840
287354456	LAT 4MM DOWN-BITING ATP CURETTE	00190776152871
287354471	LAT 4MM PUSH-PULL ATP CURETTE, CRAN	00190776152499
287354472	LAT 4MM PUSH-PULL ATP CURETTE, CAUD	00190776152505
287354481	LATERAL 4MM STRAIGHT PUSH-PULL CURETTE	00190776152192
287354484	LATERAL 4MM ANGLED PUSH-PULL CURETTE, RIGHT	00190776152338
287354485	LATERAL 4MM ANGLED PUSH-PULL CURETTE, LEFT	00190776152345
287354601	LAT 6MM STRAIGHT CURETTE	00190776152178

Catalog Number	Description	UDI
287354611	LAT 6MM ANGLED CURETTE, RIGHT	00190776152741
287354612	LAT 6MM ANGLED CURETTE, LEFT	00190776152758
287354621	LAT 6M ATP CUP CURETTE, CRAN	00190776152765
287354622	LAT 6MM ATP CUP CURETTE, CAUD	00190776152772
287354631	LAT 6MM UP-BITING CURETTE	00190776152789
287354641	LAT 6MM ANGLED UP-BITING CURETTE, RIGHT	00190776152796
287354642	LAT 6MM ANGLED UP-BITING CURETTE, LEFT	00190776152802
287354651	LAT 6MM UP-BITING ATP CURETTE	00190776152826
287354653	LAT 6MM DOWN-BITING CURETTE	00190776152482
287354654	LAT 6MM ANGLED DOWN-BITING CURETTE, RIGHT	00190776152857
287354655	LAT 6MM ANGLED DOWN-BITING CURETTE, LEFT	00190776152864
287354656	LAT 6MM DOWN-BITING ATP CURETTE	00190776152888
287354861	LAT 8MM ANGLED RING CURETTE	00190776152895
287354862	LAT 8MM ANGLED RING CURETTE, FLAT	00190776152901
287361101	LAT COBB ELEVATOR 12.7MM	00190776152055
287361102	LAT COBB ELEVATOR 20MM	00190776152062
287361103	LAT COBB ELEVATOR 12.7MM ANGLED TIP	00190776152550
287361200	LAT ANGLED COBB ELEVATOR 12.7MM, RIGHT	00190776152215
287361201	LAT ANGLED COBB ELEVATOR 12.7MM, LEFT	00190776152222
287361202	LAT 20MM ANGLED COBB ELEVATOR, RIGHT	00190776152239
287361203	LAT 20MM ANGLED COBB ELEVATOR, LEFT	00190776152246
287361300	LAT COBB ELEVATOR 12° ATP 12.7MM, RIGHT	00190776152253
287361301	LAT COBB ELEVATOR 12° ATP 12.7MM, LEFT	00190776152260
287362100	LAT 10MM STRAIGHT OSTEOTOME	00190776152024
287362101	LAT 10MM ANGLED OSTEOTOME	00190776152529
287362102	LAT 10MM STRAIGHT ATP OSTEOTOME (DINGO)	00190776152536
287362150	LAT 15MM STRAIGHT OSTEOTOME	00190776152512
287362152	LAT 15MM STRAIGHT ATP OSTEOTOME (DINGO)	00190776152543
287363104	LAT 4MM x 18MM BOX CUTTER	00190776152130
287363106	LAT 6MM x 18MM BOX CUTTER	00190776152147
287363108	LAT 8MM x 18MM BOX CUTTER	00190776152154
287363114	LAT 4MM x 18MM BOX CUTTER, ANGLED	00190776152277
287363116	LAT 6MM x 18MM BOX CUTTER, ANGLED	00190776152376
287363118	LAT 8MM x 18MM BOX CUTTER, ANGLED	00190776152383
287363134	LAT 4MM x 18MM ATP BOX CUTTER	00190776152284
287363136	LAT 6MM x 18MM ATP BOX CUTTER	00190776152628
287363138	LAT 8MM x 18MM ATP BOX CUTTER	00190776152635
287363144	LAT 4MM x 18MM ATP BOX CUTTER (DINGO)	00190776152642
287363146	LAT 6MM x 18MM ATP BOX CUTTER (DINGO)	00190776152659
287363148	LAT 8MM x 18MM ATP BOX CUTTER (DINGO)	00190776152666
287363204	LAT 4MM x 23MM BOX CUTTER	00190776152567
287363206	LAT 6MM x 23MM BOX CUTTER	00190776152574
287363208	LAT 8MM x 23MM BOX CUTTER	00190776152581
287363214	LAT 4MM x 23MM BOX CUTTER, ANGLED	00190776152598
287363216	LAT 6MM x 23MM BOX CUTTER, ANGLED	00190776152604
287363218	LAT 8MM x 23MM BOX CUTTER, ANGLED	00190776152611

Depuy Synthes ATP/Lateral Discectomy Surgical Instrumentation System Instructions

INTENDED USE

- The Disc Preparation Surgical Instruments are intended to offer a comprehensive set of surgical instruments to prepare the intervertebral disc space for interbody spinal fusion.

INTENDED USER PROFILE

- Surgical procedures should be performed only by persons having adequate training and familiarity with surgical techniques.
- Consult medical literature relative to techniques, complications and hazards prior to performance of any surgical procedure. Before using the product, all instructions regarding its safety features must be read carefully.

DEVICE DESCRIPTION

- Surgical instruments comprising fixed assemblies, simple hinged instruments and simple assemblies generally composed of medical grade stainless steels, titanium, aluminum and silicone rubber.
- Instrument case and trays may consist of different materials including stainless steels, aluminum and silicone mats.
- Devices are supplied NON-STERILE and must be inspected, cleaned and sterilized before each use.
- Devices are critical and require terminal sterilization.
- Devices are not implantable.

WARNINGS

- **To prevent damage to soft tissue or vascular structures, use the instruments carefully under fluoroscopy when crossing the disc space to release the contralateral annulus.**
- Avalign recommends thorough manual and automated cleaning of medical devices prior to sterilization. Automated methods alone may not adequately clean devices.
- Devices should be reprocessed as soon as possible following use. Instruments must be cleaned separately from cases and trays.
- All cleaning agent solutions should be replaced frequently before becoming heavily soiled.
- Prior to cleaning, sterilization and use, remove all protective caps carefully. All instruments should be inspected to ensure proper function and condition. Do not use instruments if they do not perform satisfactorily.
- The sterilization methods described have been validated with the devices in predetermined placement locations per the case and tray designs. Areas intended for specific devices shall contain only those devices.
- Risk of damage – The surgical instruments are precision devices. Careful handling is important for the accurate functioning of the devices. Improper external handling can cause the devices to malfunction.
- Use caution when handling sharp instruments to avoid injury.
- Wash the instrument case and trays with an aluminum safe, neutral pH detergent to avoid faded surface colors and deterioration of anodized surfaces.
- If a device is/was used in a patient with, or suspected of having Creutzfeldt-Jakob Disease (CJD), the device cannot be reused and must be destroyed due to the inability to reprocess or sterilize to eliminate the risk of cross-contamination.

CAUTION

Federal U.S. Law restricts this device to sale, distribution, and use, by, or on order of a physician.

LIMITATIONS ON REPROCESSING

Repeated processing has minimal effect on these instruments. End of life is normally determined by wear and damage due to use.

DISCLAIMER

It is the responsibility of the reprocessor to ensure reprocessing is performed using equipment, materials and personnel in the reprocessing facility and achieves the desired result. This requires validation and routine monitoring of the process. Any deviation by the reprocessor from the instructions provided must be properly evaluated for effectiveness and potential adverse consequences.

Reprocessing Instructions

TOOLS AND ACCESSORIES

Water	Cold Tap Water (< 20°C / 68°F) Warm Water (38° - 49°C / 100° - 120°F) Hot Tap Water (> 40°C / 104°F) Deionized (DI) or Reverse Osmosis (RO) Water (ambient)
Cleaning Agents	Neutral Enzymatic Detergent pH 6.0-8.0 i.e. MetriZyme, EndoZime, Enzol
Accessories	Assorted Sizes of Brushes and/or Pipe Cleaners with Nylon Bristles Sterile Syringes or equivalent Absorbent, Low Lint Disposable Cloths or equivalent Soaking Pans
Equipment	Medical Compressed Air Ultrasonic Cleaner (Sonicator) Automated Washer

POINT-OF-USE AND CONTAINMENT

- 1) Follow health care facility point of use practices. Keep devices moist after use to prevent soil from drying and remove excess soil and debris from all lumens, surfaces, crevices, sliding mechanisms, hinged joints, and all other hard-to-clean design features.
- 2) Follow universal precautions and contain devices in closed or covered containers for transport to central supply.


MANUAL CLEANING

- 3) Disassemble all rongeurs identified for disassembly according to these instructions:

****** Careful handling of the device in the disassembled condition is critical to avoid instrument damage.******

a. Rongeur Disassembly:

The ATP / Lateral Discectomy Surgical instrument set includes versions of the Kerrison and Pituitary rongeurs that can be disassembled to aid in reprocessing cleaning. The grey clip on the handle indicates the rongeur is a version that can be disassembled for cleaning. There may be rongeurs in the set that are not designed to be disassembled during cleaning, these rongeurs do not have the grey clips.



To disassemble the rongeurs, follow the etched instructions on instrument:

1. Squeeze the rongeur handles
2. Pull grey clip out and then turn the clip clockwise.
3. The rongeur top slide will then disassemble.

Kerrison Rongeur	Pituitary Rongeur	Note
		<p>Arrows align when the grey clip is turned and the rongeur can be disassembled</p>

The rongeurs disassemble differently depending on type. See below for how each type opens – the slide will remain attached to the rongeur body the instrument opens but the slide remains attached so that it is not lost during cleaning. The open position allows for flushing and rinsing between the slide and the bottom jaw of the rongeur.

Kerrison Rongeur	Pituitary Rongeur	Note
		<p>Disassembled device</p>
		<p>Handle the instruments carefully when in the disassembled (open) state to avoid damage to the instruments.</p>
	<p>N/A</p>	<p>For Kerrison Rongeurs ONLY: Top slide can be removed completely if necessary. When rotating a small resistance or “click” will be felt when disassembly occurs</p>

	<p>N/A</p>	<p>For Kerrison Rongeurs ONLY: If top slide was completely removed; re-attach proximal end of top slide to slot on trigger mechanism, a click will be felt when it is seated in position</p>
--	------------	---

b. Rongeur Assembly:

To re-assemble the rongeurs, (rongeurs are to be re-assembled prior to sterilization).

1. If top slide was completely removed; re-attach proximal end of top slide to slot on trigger mechanism, a click will be felt when it is seated in position.
 2. Close the slide, taking caution to carefully realign the slide with the bottom jaw of the rongeur
 3. Squeeze handles to engage top slide with bottom slide
 4. Turn the grey clip counter clockwise until clip is locked into place.
 5. Squeeze the rongeur handle to ensure complete assembly and correct function of the instrument.
- 4) Rinse devices under cold running tap water for a minimum of 3 minutes while wiping off residual soil or debris. Actuate moveable mechanisms and flush all lumens, cracks and/or crevices while rinsing.
 - 5) Prepare an enzymatic cleaning solution per manufacturer's instructions including dilution/concentration, water quality and temperature. Immerse devices and soak for a minimum of 10 minutes. While in the solution, use a soft, bristle brush to remove all traces of blood and debris from the device, paying close attention to threads, crevices, seams, and any hard to reach areas.
 - a. If the device has sliding mechanisms or hinged joints, actuate the device while scrubbing to remove trapped soil.
 - b. If the device contains a lumen, use a tight-fitting nylon brush or pipe cleaner while pushing in and out with a twisting motion to facilitate removal of debris; ensure the full diameter and depth of the lumen is accessed. Flush the lumen, three times minimum, with a syringe containing a minimum solution of 60mL.
 - 6) Remove devices and rinse/agitate in cold tap water for a minimum of 3 minutes. Actuate moveable mechanisms and flush all lumens, cracks and/or crevices while rinsing.
 - 7) Prepare a neutral detergent cleaning solution per manufacturer's instructions including dilution/concentration, water quality and temperature. Immerse devices and soak for a minimum of 5 minutes. While in the solution, use a soft, bristle brush to remove all traces of blood and debris from the device, paying close attention to threads, crevices, seams, and any hard to reach areas.
 - a. If the device has sliding mechanisms or hinged joints, actuate the device while scrubbing to remove trapped soil.
 - b. If the device contains a lumen, use a tight-fitting nylon brush or pipe cleaner while pushing in and out with a twisting motion to facilitate removal of debris; ensure the full diameter and depth of the lumen is accessed. Flush the lumen, three times minimum, with a syringe containing a minimum solution of 60mL.
 - 8) Remove devices and rinse/agitate in cold tap water for a minimum of 3 minutes. Actuate moveable mechanisms and flush all lumens, cracks and/or crevices while rinsing.
 - 9) Prepare an enzymatic cleaning solution using hot water per manufacturer's recommendations in an ultrasonic unit. Sonicate the devices for a minimum of 15 minutes using a minimum frequency of 40 kHz. It is recommended to use an ultrasonic unit with flushing attachments. Devices with lumens should be flushed with cleaning solution under the surface of the solution to ensure adequate perfusion of channels.
 - 10) Remove devices and rinse/agitate in ambient DI/RO water for a minimum of 4 minutes. Actuate moveable mechanisms and flush all lumens, cracks and/or crevices while rinsing. Flush internal lumens a minimum of 3 times with RO/DI water (minimum of 15mL) using an appropriately sized syringe. If available, use flush ports for flushing.
 - 11) Dry the device using an absorbent cloth. Dry any internal areas with filtered, compressed air.

- 12) Visually inspect the device for soil including all actuating mechanisms, cracks, crevices, and lumens. If not visibly clean, repeat steps 4-12.
- 13) Submerge device in 2-3% hydrogen peroxide. The appearance of bubbles confirms the presence of hemoglobin. Repeat steps 4-13 if bubbles appear. Adequately rinse device with DI/RO water.

AUTOMATED CLEANING

Note: All devices must be manually pre-cleaned prior to any automated cleaning process, follow steps 1-8. Steps 9-13 are optional but advised.

- 14) Transfer the devices to an automatic washer/disinfector for processing per the below minimum parameters.

Phase	Time (minutes)	Temperature	Detergent Type & Concentration
Pre-wash 1	02:00	Cold Tap Water	N/A
Enzyme Wash	02:00	Hot Tap Water	Enzyme Detergent
Wash 1	02:00	63°C / 146°F	Neutral Detergent
Rinse 1	02:00	Hot Tap Water	N/A
Purified Water Rinse	02:00	146°F / 63°C	N/A
Drying	15:00	194°F / 90°C	N/A

- 15) Dry excess moisture using an absorbent cloth. Dry any internal areas with filtered, compressed air.
- 16) Visually inspect the device for soil including all actuating mechanisms, cracks, crevices and lumens. If not visibly clean, repeat steps 4-8, 14-16.
- 17) Submerge device in 2-3% hydrogen peroxide. The appearance of bubbles confirms the presence of hemoglobin. Repeat steps 4-8, 14-17 if bubbles appear. Adequately rinse device with DI/RO water.

DISINFECTION

- Devices must be terminally sterilized (See § Sterilization).
- Aalign instruments are compatible with washer/disinfector time-temperature profiles for thermal disinfection per ISO 15883.
- Load the devices in the washer-disinfector according to the manufacturer's instructions, ensuring that the devices and lumens can drain freely.
- The following automated cycles are examples of validated cycles:

Phase	Recirculation Time (min.)	Water Temperature	Water Type
Thermal Disinfection	1	>90°C (194°F)	RI/DO Water
Thermal Disinfection	5	>90°C (194°F)	RI/DO Water

INSPECTION AND FUNCTIONAL TESTING

- Visually inspect devices for damage or wear, including sharp edges. Instruments with broken, cracked, chipped or worn features, should not be used, but should be replaced immediately.
- Verify device interfaces (junctions and threads) continue to function as intended without complications.
- Check for smooth movement of hinges. Locking mechanisms should be free of nicks.
- Lubricate before autoclaving with Instra-Lube, or a steam permeable instrument lubricant.

PACKAGING

- Only FDA cleared sterilization packaging materials should be used by the end user when packaging the devices.
- The end user should consult ANSI/AAMI ST79 or ISO 17665-1 for additional information on steam sterilization.
- **Sterilization Wrap**
 - Cases may be wrapped in a standard, medical grade sterilization wrap using the AAMI double wrap method or equivalent.
- **Rigid Sterilization Container**
 - For information regarding rigid sterilization containers, please refer to appropriate instructions for use provided by the container manufacturer or contact the manufacturer directly for guidance.

STERILIZATION

Sterilize with steam. The following are minimum cycles required for steam sterilization of Aalign devices:

1. Sterilization Wraps:

Cycle Type	Temperature	Exposure Time	Pulses	Drying Time
Prevacuum	132°C (270°F)	4 minutes	4	45 minutes

Prevacuum 134°C (273°F) 3 minutes 4 45 minutes

- The operating instructions and guidelines for maximum load configuration of the sterilizer manufacturer should be followed explicitly. The sterilizer must be properly installed, maintained, and calibrated.
- Time and temperature parameters required for sterilization vary according to type of sterilizer, cycle design, and packaging material. It is critical that process parameters be validated for each facility's individual type of sterilization equipment and product load configuration.
- A facility may choose to use different steam sterilization cycles other than the cycle suggested if the facility has properly validated the cycle to ensure adequate steam penetration and contact with the devices for sterilization. Note: rigid sterilization containers cannot be used in gravity steam cycles.
- Water droplets and visible signs of moisture on sterile packaging/wrap or the tape used to secure it may compromise the sterility of the processed loads or be indicative of a sterilization process failure. Visually check outside wrap for dryness. If there are water droplets for visible moisture observed the pack or instrument tray is considered unacceptable. Repackaging and re-sterilize the packages with visible signs of moisture.

STORAGE

- After sterilization, instruments should remain in sterilization packaging and be stored in a clean, dry cabinet or storage case.
- Care should be taken when handling devices to avoid damaging the sterile barrier.

MAINTENANCE

- Attention: Apply lubricant only on the connecting elements (locking mechanism) and moving parts.
- Discard damaged, worn or non-functional devices.

WARRANTY


- All products are guaranteed to be free from defects in material and workmanship at the time of shipping.
- Avalign instruments are reusable and meet AAMI standards for sterilization. All our products are designed and manufactured to meet the highest quality standards. We cannot accept liability for failure of products which have been modified in any way from their original design.

CONTACT

- **Notice to Patient and User:** Any serious incident that has occurred in relation to the medical devices should be reported to the manufacturer and the competent authority of the EU Member State in which the user and/or patient is established.

Manufactured by:
Avalign Technologies
8727 Clinton Park Drive
Fort Wayne, IN 46825
1-877-289-1096
www.avalign.com
product.questions@avalign.com

Authorized Representative Instruments:
Instrumed GmbH
(dba Avalign German Specialty Instruments)
78 532 Tuttlingen, Unter Buchsteig 3
Germany Phone: +49 7462 200 49 0

Authorized Representative Case and Tray:
Emergo Europe,
Prinsessegracht 20, 2514 AP
The Hague, The Netherlands

Distributed in US by:
DePuy Synthes Sales Inc. (DSSI)
325 Paramount Drive
Raynham, MA 02767 USA
www.depuysynthes.com

Distributed outside US by:
Medos International Sarl
Chemin-Blanc 38
CH-2400 LeLocle
Switzerland

Symbol	Title	Symbol	Title and Translations
	Manufacturer & Date of Manufacture		Caution
	Lot Number / Batch Code		Non-Sterile
	Catalogue Number		Federal Law (USA) restricts this device to sale by or on the order of a physician
	Consult Instructions for Use		Medical Device
	Authorized Representative in the European Community		

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